

# Biofourmis Care


## Vital Signs

You will need to manually report your vital signs on the Patient Application. See the steps below to learn how.




### Step 1

Tap a task on the home screen. For this example, tap the task "Record blood pressure".




### Step 2

Tap "Record manually instead".




### Step 3

Enter your readings from your device.



### Step 4

Tap "Save". A success message appears.



# Biofourmis Care


## Vital Signs



### Viewing Vital Signs

You can view your vital sign measurements on the "Vitals" screen by tapping the "Vitals" icon on the Biofourmis Care Patient Application's footer menu. Scroll down the "Vitals" screen to view your vital sign measurements.

Home Screen



Vitals Screen

